

Syringe 500 cc ~1s				
No	Total Volume	Perbedaan	Selisih	Deviasi
	0			
1	508	508	8	1.60%
2	1024	516	16	3.20%
3	1563	539	39	7.80%
4	2118	555	55	11.00%
5	2660	542	42	8.40%
6	3192	532	32	6.40%
7	3774	582	82	16.40%
8	4293	519	19	3.80%
9	4811	518	18	3.60%
10	5290	479	-21	4.20%
11	5836	546	46	9.20%
12	6344	508	8	1.60%
13	6905	561	61	12.20%
14	7415	510	10	2.00%
15	7942	527	27	5.40%
16	8486	544	44	8.80%
17	9032	546	46	9.20%
18	9522	490	-10	2.00%
19	10054	532	32	6.40%
20	10595	541	41	8.20%
21	11134	539	39	7.80%
22	11661	527	27	5.40%
23	12168	507	7	1.40%
24	12694	526	26	5.20%
25	13163	469	-31	6.20%

Rata-Rata Deviasi:
6.30%

Kalibrasi Airflow

Syringe 200 cc ~0.75s				
No	Total Volume	Perbedaan	Selisih	Deviasi
	0			
1	212	212	12	6.00%
2	395	183	-17	8.50%
3	615	220	20	10.00%
4	801	186	-14	7.00%
5	999	198	-2	1.00%
6	1203	204	4	2.00%
7	1401	198	-2	1.00%
8	1597	196	-4	2.00%
9	1796	199	-1	0.50%
10	2019	223	23	11.50%
11	2205	186	-14	7.00%
12	2407	202	2	1.00%
13	2594	187	-13	6.50%
14	2826	232	32	16.00%
15	2998	172	-28	14.00%
16	3224	226	26	13.00%
17	3387	163	-37	18.50%
18	3574	187	-13	6.50%
19	3762	188	-12	6.00%
20	3990	228	28	14.00%
21	4165	175	-25	12.50%
22	4358	193	-7	3.50%
23	4541	183	-17	8.50%
24	4715	174	-26	13.00%
25	4922	207	7	3.50%

Rata-Rata Deviasi:
7.72%